

Datasheet

CCT8L2 purified MaxPab rabbit polyclonal antibody (D01P)

Gene Symbol: CCT8L2

Gene Alias: CESK1, MGC118839, MGC118840

Catalog Number: H00150160-D01P

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human CCT8L2 protein.

Immunogen: CCT8L2 (NP_055221.1, 1 a.a. ~ 557 a.a) full-length human protein.

Sequence:

MDSTVPSALELPQRLALNPRESPRSPEEEEPHLLSSLA
AVQTLASVIRPCYGPGRQKFLVTMKGETVCTGCATAI
LRALELEHPAAWLLREAGQTQAENSGDGTAFVLLTE
ALLEQAEQLLKAGLPRPQLREAYATATAEVLATLPSLAI
QSLGPLEDPSWALHSVMNTHLSPMDHLTKLVAHAC
WAIKELDGSFKPERVGVICALPGGTLEDSCLLPGLAISG
KLCGQMATVLSGARVALFACPFPAHPNAPATARLSS
PADLAQFSKGSQDQLEKQVQGLAAAGINVAVVLGEVD
EETLTADKYGIVVIQARSWMEIYLLSEVLDTPLLPRLLP
PQRPGKCQRVYRQELGDGLAVVFEWECTGTPALTVV
LRGATTQGLRSAEQAVYHGIDAYFQLCQDPRLIPGAG
ATEMALAKMLSDKGSRLGPGSPAFLAFWALKYLPK
TLAENAGLAVSDVMAEMSGVHQGNLLMGVGTGIIIN
VAQEGVWDTLIVKAQGFRAVAEVLQLVTVDEIVVAKK
SPTHQEIWNPDSKKTKKHPPPVETKKILGLNN

Host: Rabbit

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 150160