

Datasheet

COMMD1 MaxPab mouse polyclonal antibody (B01)

Catalog Number: H00150684-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human COMMD1 protein.

Immunogen: COMMD1 (AAH22046.1, 1 a.a. ~ 190 a.a) full-length human protein.

Sequence:

MAAGELEGGKPLSGLLNALAQDTFHGYPGITEELLRS
QLYPEVPPEEFRPFLAKMRGILKSIASADMDFNQLEAF
LTAQTKKQGGITSDQAAVSKFWKSHKTKIRESLMNQS
RWNSGLRGLSWRVDGKSQSRHSAQIHTPVAIIELELG
KYGQSESEFLCLEFDEVKVNQILKTLSEVEESISTLISQP
N

Host: Mouse

Reactivity: Human

Applications: WB-Ti, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 150684

Gene Symbol: COMMD1

Gene Alias: C2orf5, MGC27155, MURR1

Gene Summary: COMMD1 is a regulator of copper homeostasis, sodium uptake, and NF-kappa-B (see MIM 164011) signaling (de Bie et al., 2005 [PubMed 16267171]).[supplied by OMIM]