

Datasheet

SEPT10 purified MaxPab mouse polyclonal antibody (B01P)

Catalog Number: H00151011-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human SEPT10 protein.

Immunogen: SEPT10 (NP_653311.1, 1 a.a. ~ 454 a.a) full-length human protein.

Sequence:

MASSEVARHLLFQSHMATKTCMSSQGSDDDEQIKRE
NIRSLTMSGHVGFEPLDQLVNRSIQGFCFNILCVGE
TGIGKSTLIDTLFNTNFEDYESSHFCPNVKLKAQTYELQ
ESNVQLKLTIVNTVGFQDQINKEESYQPIVDYIDAQFEA
YLQEELKIKRSLFTYHDSRIHVCLYFISPTGHSLKTLDLL
TMKNLDSKVNIIPVIAKADTVSKTELQKFKIKLMSELVSN
GVQIYQFPTDDDTIAKVNAAMNGQLPFAVVGSMDEVK
VGNKMKARQYPWGVVQVENENHCDFVKLREMLICT
NMEDLREQTHTRHYELYRRCKLEEMGFTDVGPENKP
VSVQETYEAKRHEFHGERQRKEEEMKQMFVQRVKEK
EAILKEAERELQAKFEHLKRLHQEERMKLEEKRRLLLE
EIIAFSKKKATSEIFHSQSFLATGSNLRKDKDRKNSNFL

Host: Mouse

Reactivity: Human

Applications: WB-Ti, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 151011

Gene Symbol: SEPT10

Gene Alias: FLJ11619

Gene Summary: This gene encodes a member of the septin family of cytoskeletal proteins with GTPase activity. This protein localizes to the cytoplasm and nucleus and displays GTP-binding and GTPase activity. Alternate splicing results in two transcript variants encoding different isoforms. [provided by RefSeq]