

Datasheet

NSE1 monoclonal antibody (M01), clone 1C2

Catalog Number: H00151354-M01

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a full length recombinant NSE1.

Clone Name: 1C2

Immunogen: NSE1 (AAH26346, 1 a.a. ~ 292 a.a)
full-length recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

MGNQLDRITHLNYSELPTGDPSGIEKDELRVGVAYFFS
DDEEDLDERGQPDKFGVKAPPGCTPCPESPSRHHHH
LLHQLVLNETQFSAFRGQECIFSKVSGGPQGADLSVY
AVTALPALCEPGDLLELLWLQHAPEPPAPAPHWAVYV
GGGQIIHLHQGEIRQDSLYEAGAANVGRVVNSWYRYR
PLVAELVVQNACGHLGLKSEEICWTNSESFAAWCRFG
KREFKAGGEVPAGTQPPQQYYLKVHLGENKVHTAR
FHSLEDLIREKRRIDASGRLRLVQLADLVDDKE

Host: Mouse

Reactivity: Human

Applications: ELISA, IHC-P, S-ELISA, WB-Ce, WB-Re
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG1 kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 151354

Gene Symbol: FAM84A

Gene Alias: FLJ35392, NSE1, PP11517