

Datasheet

UPP2 purified MaxPab mouse polyclonal antibody (B01P)

Catalog Number: H00151531-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human UPP2 protein.

Immunogen: UPP2 (NP_775491.1, 1 a.a. ~ 317 a.a) full-length human protein.

Sequence:

MASVIPASNRSMRSDRNTYVGKRFVHVKNPYLDLMDE
DILYHLDLGTKTHNLPAMFGDVKFVCVGGSPNRMKAF
ALFMHKELGFEAEEDIKDICAGTDRYCMYKTGPVLAI
SHGMGIPISIMLHELKLLHHARCCDVTIIRIGTSGGIGI
APGTVVITDIAVDSFFKPRFEQVILDNIVTRSTELDKELS
EELFNCSKEIPNFPTLVGHTMCTYDFYEGQGRLDGAL
CSFSREKKLDYLKRAFKAGVRNIEMESTVFAAMCGLC
GLKAAVVCVTLLDRLDCDQINLPHDVLVEYQQRPELLI
SNFIRRRRLGLCD

Host: Mouse

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 151531

Gene Symbol: UPP2

Gene Alias: UDRPASE2, UP2, UPASE2