

Datasheet

SLC36A2 monoclonal antibody (M04), clone 2H3

Catalog Number: H00153201-M04

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a partial recombinant SLC36A2.

Clone Name: 2H3

Immunogen: SLC36A2 (NP_861441, 1 a.a. ~ 71 a.a)
partial recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

MSVTKSTEGPQGAVAIKLDLMSPPEAKKLENKDSTFL
DESPSESAGLKKTGKITVFQALIHLVKGNMG TG

Host: Mouse

Reactivity: Human

Applications: ELISA, S-ELISA, WB-Ce, WB-Re
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG2b Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 153201

Gene Symbol: SLC36A2

Gene Alias: FLJ16051, MGC119658, MGC119660,
PAT2, TRAMD1