

Datasheet

PPP1R2P3 MaxPab rabbit polyclonal antibody (D01)

Catalog Number: H00153743-D01

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human PPP1R2P3 protein.

Immunogen: PPP1R2P3 (AAH66922.1, 1 a.a. ~ 205 a.a) full-length human protein.

Sequence:

MAASTASHRPIKGILKNKTSTTSSMVASAEQPRRSVDE
ELSKKSQKWDEINILATYHPADKGYGLMKIDEPSPPYH
SMMGDDDEDACRDTETTEAMAPGILAKKLAEEGLEPK
YRIQEQESSGEEDSDLSPEEREKKRQFEMRRKLHYNE
GLNIKLARQLISKDLHDDDEDEEMLETADGESMNTEES
NQGSTPSDQQQNKLRSS

Host: Rabbit

Reactivity: Human

Applications: IP, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 153743

Gene Symbol: PPP1R2P3

Gene Alias: MGC87149