

Datasheet

CNKS3 MaxPab mouse polyclonal antibody (B01)

Gene Symbol: CNKS3

Gene Alias: FLJ31349, MAGI1

Catalog Number: H00154043-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human CNKS3 protein.

Immunogen: CNKS3 (NP_775786, 1 a.a. ~ 555 a.a) full-length human protein.

Sequence:

MEPVTKWSPKQVVDWTRGLDDCLQQYVHKFEREKIN
GEQLLQISHQDLEELGVTRIGHQELVLEAVDLLCALNY
GLETDNMKNLVKLRASSHNLQNYISSRRKSPAYDGN
TSRKAPNEFLTSVVELIGAALKALLAWLDRAPFTGITDFS
VTKNKIIQLCLDLTTTVQKDCFVAEMEDKVLTVVKVLN
GICDKTIRSTTDPVMSQCACLEEVHLPNIKPGEGLGMY
IKSTYDGLHVITGTTENSPADRSQKIHAGDEVIQVNQQ
TVVGWQLKNLVKKLRENPTGVVLLKKRPTGSFNFTP
APLKNLRWKPLVQTSPPPATTQSPESTMDTSLKKEK
SAILDLYIPPPAVPYSPRDENGSEFVYGGSSKCKQPLP
GPKGSESPNSFLDQESRRRRFTIADSDQLPGYSVETNI
LPTKMREKTPSYGKPRPLSMPADGNWMGIVDPFARP
RGHGRKGEDALCRYFSNERIPPIIESSPPYRFSRPT
TERHLVRGADYIRGSRCYINSDLHSSATIPFQEEGTTK
KSGSSATKSSSTEPSSLVSWFTRLKLLTH

Host: Mouse

Reactivity: Human

Applications: Det Ab, IF, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 154043