

Datasheet

CNKS3 purified MaxPab rabbit polyclonal antibody (D01P)

Gene Symbol: CNKS3

Gene Alias: FLJ31349, MAGI1

Catalog Number: H00154043-D01P

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human CNKS3 protein.

Immunogen: CNKS3 (NP_775786.2, 1 a.a. ~ 555 a.a) full-length human protein.

Sequence:

MEPVTKWSPKQVVDWTRGLDDCLQQYVHKFEREKIN
GEQLLQISHQDLEELGVTRIGHQELVLEAVDLLCALNY
GLETDNMKNLVKLKLRASSHNQNYISSRRKSPAYDGN
TSRKAPNEFLTSVVELIGAALKALLAWLDRAPFTGITDFS
VTKNKIIQLCLDLTTTVQKDCFVAEMEDKVLTVVKVLN
GICDKTIRSTTDPVMSQCACLEEVHLPNIKPGELGMY
IKSTYDGLHVITGTTENSPADRSQKIHAGDEVIQVNQQ
TVVGWQLKNLVKKLRNP TG VVLLLKKRPTGSFNFTP
APLKNLRWKPLVQTSPPPATTQSPESTMDTSLKKEK
SAILDLYIPPPAVPYSPRDENG SFVYGGSSKCKQPLP
GPKGSESPNSFLDQESRRRRFTIADSDQLPGYSVETNI
LPTKMREKTPSYGKPRPLSMPADGNWMGIVDPFARP
RGHGRKGEDALCRYFSNERIPPIIESSPPYRFSRPT
TERHLVRGADYIRGSRCYINSDLHSSATIPFQEEGKK
KSGSSATKSSSTEP SLLVSWFTRLKLLTH

Host: Rabbit

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 154043