

Datasheet

FMR1NB polyclonal antibody (A01)

Catalog Number: H00158521-A01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length recombinant FMR1NB.

Immunogen: FMR1NB (AAH34320, 1 a.a. ~ 255 a.a)
full-length recombinant protein with GST tag.

Sequence:

MSSHRRKAKGRNRRSHRAMRVAHLELATYELAATES
NPESHHPGYEAAMADRPQPGWRESLKMRVSKPFGM
LMLSIWILLFVCYYLSYYLCSGSSYFVLANGHILPNSN
AHGQSLEEDSALEALLNFFFPTTCNLRENQVAKPCNE
LQDLSESECLRHKKCFSSSGTTSFKCFAPFRDVPKQM
MQMFGGLGAISLILVCLPIYCRSLFWRSEPADDLQRQDN
RVVTGLKKQRRKRKRKSEMLQKAARGREEHGDE

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: 50 % glycerol

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 158521

Gene Symbol: FMR1NB

Gene Alias: FLJ25736, NY-SAR-35