

Datasheet

TCEAL6 purified MaxPab mouse polyclonal antibody (B01P)

Catalog Number: H00158931-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human TCEAL6 protein.

Immunogen: TCEAL6 (AAH71675.1, 1 a.a. ~ 200 a.a) full-length human protein.

Sequence:

MEKPYNKNEGNLENEGKPEDEVEPDDEGKSDEEEKP
DAEGKTECEGKRKAEGEPGDEGQLEDKGSQEKQGGK
SEGEGKPQGEGKPASQAKPEGQPRAAEKRPAGDYVP
RKAKRKTDRTDDSPKDSQEDLQERHLSSEEMMREC
GDVSRAQEELRKKQKMGGFHWMQRDVQDPFAPRGQ
RGVVRGVRGGGRGQRGLHDIPYL

Host: Mouse

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 158931

Gene Symbol: TCEAL6

Gene Alias: Tceal3