

## Datasheet

### DYX1C1 MaxPab rabbit polyclonal antibody (D01)

**Catalog Number:** H00161582-D01

**Regulation Status:** For research use only (RUO)

**Product Description:** Rabbit polyclonal antibody raised against a full-length human DYX1C1 protein.

**Immunogen:** DYX1C1 (AAH62564.1, 1 a.a. ~ 381 a.a) full-length human protein.

**Sequence:**

MPLQVSDYSWQQTKTAVFLSLPLKGVCRDTPVFC  
TE  
NYLKVNFPFLFEAFLYAPIDESSKAKIGNDTIVFTLYK  
KEAAMWETLSVTGVDKEMMQRIREKSILQAQERAKEA  
TEAKAAAKREDQKYALSVMMKIEEEERKKIEDMKENE  
RIKATKALEAWKEYQRKAEQKKIQREEKLCQKEKQIK  
EGRKKIKYKSLTRNLASRNLPKGRNSENIFTEKLKED  
SIPAPRSVGSIKINFTPRVFPTALRESQVAEEEEWLHK  
QAEARRAMNTDIAELCDLKEEEKNPEWLKDKGNKLFA  
TENYLAAINAYNLAIRLNNKMPLLYLNRVCHLKLKNLH  
KAIEDSSKEFCSLEGIECQASEPKLSHHIPSDLHVYIQM  
A

**Host:** Rabbit

**Reactivity:** Human, Mouse

**Applications:** IF, IP, WB-Ti, WB-Tr

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Storage Buffer:** No additive

**Storage Instruction:** Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 161582

**Gene Symbol:** DYX1C1

**Gene Alias:** DYX1, DYXC1, EKN1, FLJ37882, MGC70618, RD