

Datasheet

DYX1C1 purified MaxPab rabbit polyclonal antibody (D01P)

Catalog Number: H00161582-D01P

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human DYX1C1 protein.

Immunogen: DYX1C1 (AAH62564.1, 1 a.a. ~ 381 a.a) full-length human protein.

Sequence:

MPLQVSDYSWQQTKTAVFLSLPLKGVCRD TDVFCTE
NYLKVNFPFLFEAFLYAPIDESSKAKIGNDTIVFTLYK
KEAAMWETLSVTGVDKEMMQRIREKSILQAQERAKEA
TEAKAAAKREDQKYALSVMMKIEEEERKKIEDMKENE
RIKATKALEAWKEYQRKAEEQKKIQREEKLCQKEKQIK
EGRKKIKYKSLTRNLASRN LAPKGRNSENIFTEKLKED
SIPAPRSVGSIKINFTPRVFPTALRESQVAEEEEWLHK
QAEARRAMNTDIAELCDLKEEEKNPEWLKDKGNKLFA
TENYLAAINAYNLAIRLNNKMPLLYLNRAVCHLKLKNLH
KAIEDSSKEFCSLEGIECQASEPKLSHHIPSDLHVYIQM
A

Host: Rabbit

Reactivity: Human, Mouse

Applications: IF, WB-Ti, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 161582

Gene Symbol: DYX1C1

Gene Alias: DYX1, DYXC1, EKN1, FLJ37882, MGC70618, RD