

Datasheet

PHOSPHO1 monoclonal antibody (M10A), clone 1A11

Catalog Number: H00162466-M10A

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a partial recombinant PHOSPHO1.

Clone Name: 1A11

Immunogen: PHOSPHO1 (NP_848595, 168 a.a. ~ 266 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

CARCPANMCKHKVLS DYLRERAH DGVHFERLFYVGD
GANDFCPMGLLAGGDVAFPRRGYP MHR LIQE AQKAE
PSSFRASVVPWETAADVRLHLQQVLKS

Host: Mouse

Reactivity: Human, Rat

Applications: ELISA, WB-Ce, WB-Re, WB-Ti
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Isotype: IgG1 Kappa

Storage Buffer: In ascites fluid

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 162466

Gene Symbol: PHOSPHO1

Gene Alias: -