

## Datasheet

### OXER1 purified MaxPab mouse polyclonal antibody (B01P)

**Catalog Number:** H00165140-B01P

**Regulation Status:** For research use only (RUO)

**Product Description:** Mouse polyclonal antibody raised against a full-length human OXER1 protein.

**Immunogen:** OXER1 (NP\_683765.1, 1 a.a. ~ 423 a.a) full-length human protein.

**Sequence:**

MLCHRGGQLIVPIIPLCPEHSCRGRRLLQNLLSGPWPK  
QPMELHNLSSPSPSLSSSVLPPSFSPSPSSAPSAFTTV  
GGSSGGPCHPTSSSLVSAFLAPILALEFVLGLVGNSLA  
LFIFCIHTRPWTSTNTVFLVSLVAADFLLISNLPLRVDYYL  
LHETWRFGAAACKVNLFMLSTNRTASVVFLTAIALNRY  
LKVVQPHHVLSRASVGAAARVAGGLWVGILLNGHLL  
LSTFSGPSCLSYRVGKPSASLRWHQALYLEFFLPLA  
LILFAIVSIGLTIRNRGLGGQAGPQRAMRVLAMVVAVYT  
ICFLPSIIFGMASMAFVWLSACRSLDLCTQLFHGSLAFT  
YLNVLDPVLYCFSSPNFLHQSRALLGLTRGRQGPVS  
DESSYQPSRQWRYREASRKAEAGKLVQGEVSLEKE  
GSSQG

**Host:** Mouse

**Reactivity:** Human

**Applications:** WB-Ti, WB-Tr

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Storage Buffer:** In 1x PBS, pH 7.4

**Storage Instruction:** Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 165140

**Gene Symbol:** OXER1

**Gene Alias:** GPCR, GPR170, TG1019