

Datasheet

ZNF509 purified MaxPab mouse polyclonal antibody (B01P)

Catalog Number: H00166793-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human ZNF509 protein.

Immunogen: ZNF509 (AAH16477.1, 1 a.a. ~ 248 a.a) full-length human protein.

Sequence:

MQRKLVKHIRHTGERPYSCSACGKCFGSGDLRRH
VRTHTGEKPYTCEICNKCFTSAVLRRHKMHCKAGD
ESPDVLEELSQAIENTSDLEKSQSSDSFSQDTSVTLMPV
SVKLPVHPVENSVEFDSHSGGSYCKLRSMIQPHGVS
DQEKLSLDPGKLAKPQMQQTQPQAYAYSVDVTPAGG
EPLQADGMAMIRSSLAALDNHGGDPLGSRASSTTYRN
SEGQFFSSMTLWGLAMKTLQENELDQ

Host: Mouse

Reactivity: Human

Applications: IF, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 166793

Gene Symbol: ZNF509

Gene Alias: FLJ38559, MGC126279, MGC126280