

Datasheet

DCP2 purified MaxPab mouse polyclonal antibody (B01P)

Catalog Number: H00167227-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human DCP2 protein.

Immunogen: DCP2 (AAH64593.1, 1 a.a. ~ 385 a.a) full-length human protein.

Sequence:

METKRVEIPGSVLDDLCR FILHIPSEERDN AIRVCFQI
ELAHWFYLD FYMQNTPGLPQCGIRDFAKAVFSHC PFL
LPQGEDVEKVLDEWKEYKMGVPTYGAIIDETLE NVLL
VQGYLAKSGWGFPGKGVNKEEAPHDCAAREVFEETG
FDIKDYICKDDYIELRINDQLARLYIIPGIPKDTKFNPKTR
REIRNIEWFSIEKLPCHRNDMTPKSKLGLAPNKFFMAIP
FIRPLRDWLSRRFGDSSDSDNGFSSTGSTPAKPTVEK
LSRTKFRHSQQLFDPDGSPGDQWVKHRQPLQQKPYNN
HSEMSDLLKGKKCEKKLHPRKLQDNFETDAVYDLPSS
SEDQLEHAEGQPVACNGHCKFPFSSRAFLSFKFDHN
AIMKILDL

Host: Mouse

Reactivity: Human

Applications: IF, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 167227

Gene Symbol: DCP2

Gene Alias: FLJ33245, NUDT20

Gene Summary: DCP2 is a key component of an

mRNA-decapping complex required for removal of the 5-prime cap from mRNA prior to its degradation from the 5-prime end (Fenger-Gron et al., 2005 [PubMed 16364915]).[supplied by OMIM]