

Datasheet

PRSS35 purified MaxPab mouse polyclonal antibody (B01P)

Catalog Number: H00167681-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human PRSS35 protein.

Immunogen: PRSS35 (NP_699193, 1 a.a. ~ 413 a.a) full-length human protein.

Sequence:

MENMLLWLIFFTPGWTLIDGSEMEWDFMWHLRKVPRI
VSERTFHLTSPAFAEADAKMMVNTVCGIECQKELPTPSL
SELEDYLSYETVFENGTRTLTRVKVQDLVLEPTQNITT
KGVSVRRKRQVYGTDSRFSILDKRFLTNPFPSTAVKLS
TGCSGILISPQHVLTAACHVHDGKDYVKGSKKLRVGLL
KMRNKSGGKKRRGSKRSRREASGGDQREGTREHLQ
ERAKGGRRRKSGRGQRIAEGRPSFQWTRVKNTHIP
KGWARGGMGDATLDYDYALLELKRAHKKKYMELGISP
TIKKMPGGMIHFSGFDNDRADQLVYRFCSVDESNDL
LYQYCDAESGSTGSGVYLRKDPDKKNWKRKIIAVYS
GHQWVDVHGVQKDYNVAVRITPLKYAQICLWIHGND
NCAYG

Host: Mouse

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 167681

Gene Symbol: PRSS35

Gene Alias: C6orf158, MGC46520, dJ223E3.1