

Datasheet

ZNF92 purified MaxPab rabbit polyclonal antibody (D01P)

Gene Alias: HPF12, HTF12, TF12

Catalog Number: H00168374-D01P

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human ZNF92 protein.

Immunogen: ZNF92 (AAH36439.1, 1 a.a. ~ 517 a.a) full-length human protein.

Sequence:

MVDKTPVMCSHFAQDVWPEHSIKDSFQKVILRTYGKY
GHENLQLRKDHKSVDACKVYKGGYNGLNQCLTTTDS
KIFQCDKYVKVFHKFPNVNRNKIRHTGKKPFCKNKG
KSFCMLSQLTQHKKIHTREYSYKCEECGKAFNWSSTL
TKHKIIHTGEKPYKCEECGKAFNRSSNLTKHKIIHTGEK
PYKCEECGKAFNRSSNLTKHKRIHTEKPYKCEECGK
AFNQFSILNKHKRIHMEDKPYKCEECGKAFRVFSILKK
HKIIHTGEKPYKCEECGKAFNQFSNLTKHKIIHTGEKPY
KCDECGKAFNQSSNLTKHKRIHTGEKPYKCEECGKAF
KQSSLTTEHKIIHTGEKPYKCEKCGKAFSWSSAFTKHK
RNHMEDKPYKCEECGKAFSVFSTLTKHKIIHTREKPYK
CEECGKAFNRSSIFTKHKIIHTEGKSYKCEKCGNAFNQ
SSNLTKARKIIYTGEKPYKYEECDKAFNKFSTLITHQIIYT
GEKPCKHECGRAFKNSSNYTKEKLQT

Host: Rabbit

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 168374

Gene Symbol: ZNF92