

Datasheet

ZNF169 MaxPab mouse polyclonal antibody (B01)

Entrez GeneID: 169841**Gene Symbol:** ZNF169**Gene Alias:** MGC51961**Catalog Number:** H00169841-B01**Regulation Status:** For research use only (RUO)**Product Description:** Mouse polyclonal antibody raised against a full-length human ZNF169 protein.**Immunogen:** ZNF169 (AAH47702.1, 1 a.a. ~ 590 a.a) full-length human protein.**Sequence:**

MAFRDVAVAFQTQKEWKLLSSAQRTLYREVMLENYSHL
VSLGIAFSKPKLIEQLEQGDEPWRENEHLLDLCPEPR
TEFQPSFPHLVAFSSSQLLRQYALSGHPTQIFPSSSAG
GDFQLEAPRCSSEKGESGETEGPDSSLRKRPSRISRT
FFSPHQGDPVEWVEGNREGGTDLRLAQRMSSLGSD
TMLKGADTSESGAVIRGNRYRLGLSKKSSLSFQKHHV
CPECGRGFCQQRSDLIKHQRTHTEKPYLCPECGRRF
SQKASLSIHQRKHSGEKPYVCRECGRHFRYTSSLTNH
KRIHSGERPFVQCQECGRGFRQKIALLLHQRTHLEEKPF
VCPECGRGFCQKASLLQHQSSTGERPFLCLEGRS
FRQQSLLLSHQVTHSGEKPYVCAECGHSFRQKVTLIR
HQRTHTEKPYLCPCQGRGFSQKVTLIGHQRTHTE
KPYLCPCDGRGFGQKVTLIRHQRTHTEKPYLCPCG
RAFGFKSLLTRHQRTHSEELYVDRVCGQLGQKSHL
ISDQRTHSGEKPCICDECGRGFGFKSALIRHQRTHSG
EKPYVCRECGRGFSQKSHLHRHRTKSGHQLLPQEV
F

Host: Mouse**Reactivity:** Human**Applications:** IF, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols**Storage Buffer:** No additive**Storage Instruction:** Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.