

Datasheet

NUDT10 purified MaxPab rabbit polyclonal antibody (D01P)

Catalog Number: H00170685-D01P

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human NUDT10 protein.

Immunogen: NUDT10 (NP_694853.1, 1 a.a. ~ 164 a.a) full-length human protein.

Sequence:

MKCKPNQTRTYDPEGFKKRAACLCFRSEREDEVLLVS
SSRYPDRWIVPGGGMEPEEEPGGAHVREYEEAGVK
GKLGRLLGVFEQNQDPKHRTYVYVLTVELLEDWEDS
VSGRKREWFKVEDAIKVLQCHKPVHAEYLEKLKLGGS
PTNGNSMAPSSPDSDP

Host: Rabbit

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 170685

Gene Symbol: NUDT10

Gene Alias: APS2, DIP3a, hDIP3alpha

Gene Summary: NUDT10 belongs to a subgroup of phosphohydrolases that preferentially attack diphosphoinositol polyphosphates (Hidaka et al., 2002 [PubMed 12105228]).[supplied by OMIM]