

Datasheet

CCDC131 purified MaxPab mouse polyclonal antibody (B01P)

Catalog Number: H00196441-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human CCDC131 protein.

Immunogen: CCDC131 (ENSP00000307825, 1 a.a. ~ 358 a.a) full-length human protein.

Sequence:

MATADTPAPASSGLSPKEEGELEDGEISDDDNNSQIR
SRSSSSSSGGGLLPYPRRRPPHSARGGGSGGGGGS
SSSSSSSQQLRNFSRSRHASERGHRLGPSSYRPKE
PFRSHPPSVRMPSSSLSESSPRPSFWERSHLALDRFR
FRGRPYRGGSRSRGRGVGERGGKPGCRPPLGGG
AGSGFSSSQSWREPSPPRKSSKSFGFRSPSRKQNYSS
KNENCVEETFEDLLLKYKQIQLELECINKDEKLALSSKE
ENVQEDPKTLNFEDQTSTDNVSITKDSSKEVAPEEKT
QVKTFQAFELKPLRQKLTLPGDKNRLKKVKDGAKPLSL
KSDTTDSSQGIPYRVKEGFTPIPLKFSFA

Host: Mouse

Reactivity: Human

Applications: IF, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 196441

Gene Symbol: ZFC3H1

Gene Alias: CCDC131, DKFZp686A0722, KIAA0546, MGC23401, MGC90200, PSRC2