

Datasheet

LCTL purified MaxPab mouse polyclonal antibody (B01P)

Gene Symbol: LCTL

Gene Alias: FLJ33279, KLG, KLPH

Catalog Number: H00197021-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human LCTL protein.

Immunogen: LCTL (AAI52991.1, 1 a.a. ~ 567 a.a)
full-length human protein.

Sequence:

MKPVVWVATLLWMLLLVPRLGAARKGSPEEASFYYGTF
PLGFSWGVGSSAYQTEGAWDQDGKGPSIWDVFTHS
GKGKVLGNETADVACDGYKQVEDIILLRELHVNHYRF
SLSWPRLPTGIRAEQVNKKGIEFYSDLIDALLSSNITPI
VTLHHWDLPQLLQVKYGGWQNVSMANYFRDYANLCF
EAFGDRVKHWTFSDPRAMAEKGYETGHHAPGLKLR
GTGLYKAAHHIIKAHAKAWHSYNTTWRSKQQGLVGIS
LNCDWGEPVDISNPKDLEAAERYLQFCLGWFIANPIYA
GDYPQVMKDYIGRKSAREEQGLEMSRLPVFSLQEKSIIK
GTSDFLGLGHFTTRYITERNYPSRQGPPSYQNDRLIEL
VDPNWPDLGSKWLYSVPWGFRRLLNFAQTQYGDPPPI
YVMENGASQKFHCTQLCDEWRIQYLKGYINEMLKAIK
DGANIKGYTSWLLDKFEWEKGYSDRYGFYYVEFND
RNKPRYPKASVQYYKKIIIANGFPNPREVESWYLKALE
TCSINNQMMLAAEPLLSHMQMVTIIVPTVCSLCVLITAV
LLMLLLRRQS

Host: Mouse

Reactivity: Human

Applications: Det Ab, WB-Ti, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 197021