

Datasheet

ZNF553 MaxPab mouse polyclonal antibody (B03)

Catalog Number: H00197407-B03

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human ZNF553 protein.

Immunogen: ZNF553 (NP_689865.1, 1 a.a. ~ 618 a.a) full-length human protein.

Sequence:

MERAVEPWGPD LHRPEERE PQRGARTGLGSENVISQ
PNEFEHTPQEDDLGFKEEDLAPDHEVGNVSLKPEGIQ
NWDDLWVQREGLGKPQPRDRGPRLLGEPRWQASS
DRAAVCGECGKSFRQMSDLVKHQRTHTGEKPYKCGV
CGKGFGDSSARIKHQRTHSGEKPYRARPPAQGPPKIP
RSRIPAGERPTICGECGKSFRQSSDLVKHQRTHTGEK
PYKCGICGKGFGDSSARIKHQRTHRGEQPPRPVVPRR
QPSRAATAATQGPKAQDKPYICTDCGKRFLVSCSLLS
HQRSHLGPKPFGCDVCGKEFARGSDLVKHLRVHTGE
KPYLCPECGKGFADSSARVKHLRTHSGERPHACPEC
DRTFSLSSTLLRHRLTHMEPQDFSFGYPLPALIPSP
PPPLGTSPPLTPRSPSHSGEPFGLPGLEPEPGGPQAG
EPPPLAGDKPHKCPECGKGFRRSSDLVKHHRVHTG
EKPYLCPECGKGFADSSARVKHLRTHRGERARPPPP
STLLRPHNPPGPVPMAPRPRVRAQPSGPSQPHVCGF
CGKEFPRSSDLVKHRRHTHTGEKPYKCAECGKFGDS
SARIKHQRGHLVLT PFGIGDGRARPLKQEAATGLE

Host: Mouse

Reactivity: Human

Applications: Det Ab, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 197407

Gene Symbol: ZNF48

Gene Alias: DKFZp762K013, FLJ31751, MGC43952, ZNF553