

Datasheet

MGC27277 monoclonal antibody (M01), clone 3C12

Catalog Number: H00200095-M01

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a full length recombinant MGC27277.

Clone Name: 3C12

Immunogen: MGC27277 (AAH42869, 1 a.a. ~ 189 a.a)
full-length recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

MDKEETKTKPRLLRYEEKKYEDVKPLETQPAEIAEKET
LEYKTVRTFSESLKSEKTEDYLRESIIQQHMQVSEPEASL
KEGKGSRRKKDQTHACPNVRKARPVSYDRTEPKDDD
VIRNIIRLREKLGWQTILPQHSLKYGSSKIAIQKITLKKPL
EDDGEFVYCLPRKSPKSLNPNYDLQVINIVGSVKE

Host: Mouse

Reactivity: Human

Applications: ELISA, S-ELISA, WB-Re
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG1 kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 200095

Gene Symbol: C1orf67

Gene Alias: MGC149665, MGC149666, MGC27277,
MGC51214