

Datasheet

RTKN2 monoclonal antibody (M03), clone 2C2

Catalog Number: H00219790-M03

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a full-length recombinant RTKN2.

Clone Name: 2C2

Immunogen: RTKN2 (AAH25765, 1 a.a. ~ 163 a.a)
full-length recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

MEGPSLRG PALRLAGLPTQQDCNIQEKIDLEIRMREGI
WKLLSLSTQKDQVLHAVKNLMVCNARLMAYTSELQKL
EEQIANQTGRCDVKFESKERTACKGKIAISDIRIPLMWK
DSDHFSNKERSRRYAIFCLFKMGANVFDTDVVNVDKTI
TDICFENVIL

Host: Mouse

Reactivity: Human

Applications: ELISA, IHC-P, S-ELISA, WB-Re
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG2a Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 219790

Gene Symbol: RTKN2

Gene Alias: DKFZp686J10120, PLEKHK1, bA531F24.1