

Datasheet

MRGPRF MaxPab mouse polyclonal antibody (B01)

Catalog Number: H00219928-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human MRGPRF protein.

Immunogen: MRGPRF (AAH16964, 1 a.a. ~ 343 a.a)
full-length human protein.

Sequence:

MAGNCSWEAHPGNRNKMCPGLSEAPELYSRGFLTIE
QIAMLPPPAVMNYIFLLLCLGLVGNGLVLWFFGFSIKR
NPFISYFLHLASADVGLFSKAVFSILNTGGFLGTFADY
IRSVCRVLGLCMFLTGVSLPAVSAERCASVIFPAWYW
RRRPKRLSAVVCALLWVLSLLVTCLHNYFCVFLGRGA
PGAACRHMDIFLGILLFLLCCPLMVLPCALILHVECRA
RRRQRS AKLNHVILAMVSVFLVSSIYLGIDWFLFWVFQI
PAPFPEYVTDLCICINSSAKPIVYFLAGRDKSQRLWEPL
RVVFQRALRDGAELGEAGGSTPNTVTMEMQCPPGNA
S

Host: Mouse

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 219928

Gene Symbol: MRGPRF

Gene Alias: GPR140, GPR168, MGC21621, RTA, mrgF