

Datasheet

LOC220594 polyclonal antibody (A01)

Catalog Number: H00220594-A01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a partial recombinant LOC220594.

Immunogen: LOC220594 (NP_665808, 255 a.a. ~ 354 a.a) partial recombinant protein with GST tag.

Sequence:

HIKPIYNLYAISCHSGILSGGHVVTYAKNPCKWYCYN
GSICEEHPDEIDTDSAYILFYEQQRIDYAQFLPKIDGK
KMADTSSMDEDFESDYEKYCVLQ

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Ce, WB-Re

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: 50 % glycerol

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 220594

Gene Symbol: LOC220594

Gene Alias: FLJ98241