

Datasheet

CD5L purified MaxPab rabbit polyclonal antibody (D01P)

Catalog Number: H00000922-D01P

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human CD5L protein.

Immunogen: CD5L (NP_005885.1, 1 a.a. ~ 347 a.a) full-length human protein.

Sequence:

MALLFSLILAICTRPGFLASPSGVRLVGGLHRCEGRVE
VEQKGQWGTVCDDGWDIKDVAVLCRELGCGAASGTP
SGILYEPPAEKEQKVLIQSVSCTGTEDTLAQCEQEEVY
DCSHDEDAGASCENPESSFSPVPEGVRLADGPGHCK
GRVEVKHQNQWYTVCTGWSLRAAKVVCRQLGCGR
AVLTQKRCNKHAYGRKPIWLSQMSCSGREATLQDCP
SGPWGKNTCNHDEDTWVECEDPFDLRLVGGDNLCS
GRLEVLHKGWGSVCDDNWGEKEDQVVCKQLGCGK
SLSPSFRDRKCYGPGVGRIWLDNVRCSGEEQSLEQC
QHRFWGFHDCTHQEDVAVICSG

Host: Rabbit

Reactivity: Human, Mouse

Applications: WB-Ce, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 922

Gene Symbol: CD5L

Gene Alias: AIM, API6, PRO229, SP-ALPHA, Spalpha