

Datasheet

SIGLEC6 monoclonal antibody (M02), clone 2G6

Catalog Number: H00000946-M02

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a partial recombinant SIGLEC6.

Clone Name: 2G6

Immunogen: SIGLEC6 (NP_001236, 371 a.a. ~ 453 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

KTRRKAAQPVQNTDDVNPVMVSGSRGHQHQFQTGI
VSDHPAEAGPISEDEQELHYAVLHFHKVQPQEPKVTD
TEYSEIKIHK

Host: Mouse

Reactivity: Human

Applications: ELISA, IP, WB-Re, WB-Tr
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Isotype: IgG2b Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 946

Gene Symbol: SIGLEC6

Gene Alias: CD327, CD33L, CD33L1, CDw327, OBPP1, SIGLEC-6