

Datasheet

BBX purified MaxPab mouse polyclonal antibody (B01P)

Catalog Number: H00056987-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human BBX protein.

Immunogen: BBX (NP_064620.2, 1 a.a. ~ 911 a.a) full-length human protein.

Sequence:

MKGSNRNKDHS AEGEGVGKRPKRKCLQWHPLLAKKL
LDFSEEEEEDEEEDIDKVQLLGADGLEQDVGETEDD
ESPEQRARRPMNAFLLFCKRHRSLVRQEHPRLDNRG
ATKILADWWAVLDPKEKQKYTDMAKEYKDAFMKANP
GYKWCPTTNKPKVSPTPTVNPRKKLWAFPSDSSRD
PSPKKAKEEMPQLNFGMADPTQMGGLSMLLLAGEH
ALGTPEVSSGTCRPDVSESPELRQKSPLFQFAEISSST
SHSDASTKQCQTSALFQFAEISSNTSQLGGAEPVKRC
GKSALFQLAEMCLASEGMKMEESKLIKAKESDGGRIK
ELEKGKEEKEIKMEKTDETRLQKEAEFEKSAKENLRDS
KELRNFEALQIDDIMAIKMEDPKEIRKEELEEDHKCSHF
PDFSYSASSKIIISDVPSRKDHMCHPHGIMIIEDPAALNK
PEKLLKKKKKSKMDRHGNDKSTPKKTCKKRQSSESDI
ESVIYTIEAVAKGDWGIEKLGDTPRKKVVRTSSSGKGSIL
DAKPPKKKVKVSREKKMSKEKSSDTTKESRPPDFISISA
SKNISGETPEGIKAEP LTPMEDALPPSLSGQAKPEDSD
CHRIETCGSRKRSERCKGALYKTLVSEGMLTSLRAN
VDRGKRSSGKGNSSDHEGCWNEESWTFSQSGTSGS
KKFKKTKPKEDCLLGS AKLDEEFKFKFNSLPQYSPVTF
DRKCVPVPRKKKKTGNVSSEPTKTSKSGDKWSNKQ
LFLDAIHPTEAIFSED RNTMEPVHKVKNIPSIFNTPEPTT
TQEPLVGSQKRKARKTKITHLVRTADGRVSPAGGTLD
DKPKEQLQRSLPKATETDCNDKCSHNTEVGETRSSTP
EMPAVS AFFSLAALAEVAAMENVHRGQRSTPLTHDG
QPKEMPQAPVLISCADQ

Host: Mouse

Reactivity: Human

Applications: Det Ab, WB-Ti, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 56987

Gene Symbol: BBX

Gene Alias: HBP2, HSPC339, MDS001