

Datasheet

KIF15 polyclonal antibody (A01)

Catalog Number: H00056992-A01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a partial recombinant KIF15.

Immunogen: KIF15 (NP_064627, 1279 a.a. ~ 1385 a.a) partial recombinant protein with GST tag.

Sequence:

NKEMECLRMTDEVERTQTLESKAFQEKEQLRSKLEEM
YEERERTSQEMEMLRKQVECLAEENGKLVGHQNLHQ
KIQYVRLKKENVRLAEETEKLRANVFLKEKKR

Host: Mouse

Reactivity: Human

Applications: ELISA

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: 50 % glycerol

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 56992

Gene Symbol: KIF15

Gene Alias: FLJ25667, HKLP2, KNSL7, NY-BR-62