

Datasheet

SAS10 MaxPab mouse polyclonal antibody (B01)

Catalog Number: H00057050-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human SAS10 protein.

Immunogen: SAS10 (NP_065101.1, 1 a.a. ~ 479 a.a) full-length human protein.

Sequence:

MVGRSRRRGAAKWA AVRAKAGPTLTDENGDDLGLPP
SPGDTSY YQDQVDDFHEARSRAALAKGWNEVQSGD
EEDGEEEEEEVLALDMDDEDEDGGNAGEEEEEENA
DDDGGSSVQSEAEASVDP SLWGQRKKLYD TDYGS
KSRGRQSQQEAE EEEEEEEEEEAQIIQRRLAQALQEDD
FGVAWVEAFAPVPQVDEAETRVVKDLAKVSVKEK LK
MLRKESPELLELIEDLKVKLTEVKDELEPLLELVEQGIIP
PGKGSQYLRTKYNL YLNYCSNISFYLIKARRVPAHGH
PVIERLV TYRNLINKLSVVDQKLSSEIRHLLTLKDDAVK
KELIPKAKSTKPKPKSVSKTSAAACAVTDLSDDSDFDE
KAKLKYYKEIEDRQKLKRKKEENSTEEQAL EDQNAKR
AITYQIAKNRGLTPRRKKIDRNPRVKHREKFRRAKIRR
RGQVREVRKEEQRYSGELSGIRAGVKKSIKLK

Host: Mouse

Reactivity: Human

Applications: IF, WB-Ce, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 57050

Gene Symbol: UTP3

Gene Alias: CRL1, CRLZ1, FLJ23256, SAS10