

Datasheet

AGTRAP polyclonal antibody (A01)

Catalog Number: H00057085-A01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a partial recombinant AGTRAP.

Immunogen: AGTRAP (NP_065083, 108 a.a. ~ 159 a.a) partial recombinant protein with GST tag.

Sequence:

HYRERGGELLVHTGFLGSSQDRSAYQTIDSAEAPAD
PFAVPEGRSQDARGS

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: 50 % glycerol

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 57085

Gene Symbol: AGTRAP

Gene Alias: ATRAP, MGC29646

Gene Summary: This gene encodes a transmembrane protein localized to the plasma membrane and perinuclear vesicular structures. The gene product interacts with the angiotensin II type I receptor and negatively regulates angiotensin II signaling. Alternative splicing of this gene generates multiple transcript variants encoding different isoforms. [provided by RefSeq]