

Datasheet

RCAN3 polyclonal antibody

Catalog Number: PAB4905

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of RCAN3.

Immunogen: A synthetic peptide (conjugated with KLH) corresponding to N-terminus of human RCAN3.

Host: Rabbit

Reactivity: Human, Mouse

Applications: ELISA, WB-Ti

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Purification: Protein G purification

Recommend Usage: ELISA (1:1000)

Western Blot (1:100-500)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.09% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -20°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 11123

Gene Symbol: RCAN3

Gene Alias: DSCR1L2, MCIP3, RCN3, hRCN3