

Datasheet

Myeloid Specific Antigen monoclonal antibody, clone BM-2

Catalog Number: MAB5183

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against native Myeloid Specific Antigen.

Clone Name: BM-2

Immunogen: Native purified from pokeweed mitogen
stimulated human peripheral blood lymphocyte nuclei.

Host: Mouse

Reactivity: Human

Applications: IHC-Fr, IHC-P

(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Specificity: BM-2 recognizes a cytoplasmic antigen
expressed in mature human granulocytes (polys)
residing in lymphoid and non-lymphoid tissues. It does
not react with any other cell type in human tissues.

Form: Liquid

Isotype: IgG1

Recommend Usage: The optimal working dilution
should be determined by the end user.

Storage Buffer: In PBS (0.05% sodium azide)

Storage Instruction: Store at 4°C.