

Datasheet

IRS1 (phospho S312) polyclonal antibody

Catalog Number: PAB5738

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic phosphopeptide of IRS1.

Immunogen: Synthetic phosphopeptide corresponding to residues surrounding S312 of human IRS1.

Host: Rabbit

Reactivity: Human, Mouse, Rat

Applications: WB-Ce

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Recommend Usage: Western Blot (1:500-1:1000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, 150 mM NaCl, pH 7.4 (50% glycerol, 0.02% sodium azide)

Storage Instruction: Store at -20°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 3667

Gene Symbol: IRS1

Gene Alias: HIRS-1