

Datasheet

GSK3A monoclonal antibody, clone SBGSK3a4

Catalog Number: MAB5634

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against recombinant GSK3A.

Clone Name: SBGSK3a4

Immunogen: Recombinant protein corresponding to
GSK3A.

Host: Mouse

Reactivity: Human, Mouse, Rat

Applications: IP, WB

(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Specificity: human, Mouse, and Rat GSK-3a.

Form: Liquid

Isotype: IgG2a

Recommend Usage: The optimal working dilution
should be determined by the end user.

Storage Buffer: In 100 mM BBS, pH 8.2

Storage Instruction: Store at 4°C.