

Datasheet

Phosphoserine polyclonal antibody

Catalog Number: PAB13796

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against Phosphoserine.

Immunogen: Phosphoserine conjugated with KLH and phosvitin mixture.

Host: Rabbit

Applications: ELISA, IP, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: Recognizes proteins phosphorylated on serine residues. Does not cross-react with phosphotyrosine.

Form: Liquid

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: In TBS (0.1% sodium azide)

Storage Instruction: Store at -20°C.

Aliquot to avoid repeated freezing and thawing.