

Datasheet

LAMP2 polyclonal antibody

Catalog Number: PAB13805

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of LAMP2.

Immunogen: A synthetic peptide corresponding to internal region 17 amino acids of human LAMP2.

Host: Rabbit

Reactivity: Human, Mouse

Applications: WB-Ce

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: Recognizes LAMP-2.

Form: Liquid

Recommend Usage: Western Blot (1-2 ug/mL)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.02% sodium azide)

Storage Instruction: Store at 4°C.

Entrez GeneID: 3920

Gene Symbol: LAMP2

Gene Alias: CD107b, LAMPB, LGP110

Gene Summary: The protein encoded by this gene is a member of a family of membrane glycoproteins. This glycoprotein provides selectins with carbohydrate ligands. It may play a role in tumor cell metastasis. It may also function in the protection, maintenance, and adhesion of the lysosome. Alternative splicing of this gene results in multiple transcript variants encoding distinct proteins. [provided by RefSeq]