

## Datasheet

### TP53INP1 polyclonal antibody

**Catalog Number:** PAB13806

**Regulation Status:** For research use only (RUO)

**Product Description:** Rabbit polyclonal antibody raised against synthetic peptide of TP53INP1.

**Immunogen:** A synthetic peptide corresponding to N-terminus 14 amino acids of human TP53INP1.

**Host:** Rabbit

**Reactivity:** Human, Mouse, Rat

**Applications:** IHC-P, WB-Ti  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at  
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Specificity:** Recognizes p53DINP1.

**Form:** Liquid

**Recommend Usage:** Western Blot (1-2 ug/mL)  
Immunohistochemistry (2 ug/mL)  
The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS (0.02% sodium azide)

**Storage Instruction:** Store at 4°C.

**Entrez GeneID:** 94241

**Gene Symbol:** TP53INP1

**Gene Alias:** DKFZp434M1317, FLJ22139, SIP,  
TP53DINP1, TP53INP1A, TP53INP1B, Teap, p53DINP1