

Datasheet

Rubella Virus Virions polyclonal antibody

Catalog Number: PAB13818

Regulation Status: For research use only (RUO)

Product Description: Goat polyclonal antibody raised against Rubella Virus Virions.

Immunogen: Rubella Virus, strain HPV77.

Host: Goat

Reactivity: Viruses

Applications: HI, IF

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: Purified virions. Negative against VERO cells by direct immunofluorescence.

Form: Liquid

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: In 10 mM PBS, pH 7.2 (0.1% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -20°C.

Aliquot to avoid repeated freezing and thawing.