

## Datasheet

### ***Staphylococcus aureus* polyclonal antibody**

**Catalog Number:** PAB13827

**Regulation Status:** For research use only (RUO)

**Product Description:** Rabbit polyclonal antibody raised against *Staphylococcus aureus*.

**Immunogen:** Native from *Staphylococcus aureus* (ATCC #27660).

**Host:** Rabbit

**Reactivity:** Staphylococcus

**Applications:** IF

(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Specificity:** Soluble and structural antigens of the whole bacterium. Antiserum is not absorbed for related microorganisms.

**Form:** Liquid

**Recommend Usage:** The optimal working dilution should be determined by the end user.

**Storage Buffer:** In 10 mM PBS, pH 7.2 (0.1% sodium azide)

**Storage Instruction:** Store at 4°C. For long term storage store at -20°C.  
Aliquot to avoid repeated freezing and thawing.