

Datasheet

***Staphylococcus aureus* Enterotoxin B (SEB) polyclonal antibody**

Catalog Number: PAB13829

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against *Staphylococcus aureus* Enterotoxin B.

Immunogen: Purified *Staphylococcus aureus* Enterotoxin B.

Host: Rabbit

Reactivity: Staphylococcus

Applications: ID

(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: Staphylococcus aureus Enterotoxin B.
Non-reactive with other Staph. toxins.

Form: Liquid

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: In 10 mM PBS, pH 7.2 (0.1% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -20°C.
Aliquot to avoid repeated freezing and thawing.