

Datasheet

***Streptococcus* Group A polyclonal antibody**

Catalog Number: PAB13830

Regulation Status: For research use only (RUO)

Product Description: Goat polyclonal antibody raised against *Streptococcus* Group A.

Immunogen: *Streptococci* Group A.

Host: Goat

Reactivity: Staphylococcus

Applications: Cap Ab, Det Ab, ELISA, LFIA, S-ELISA
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: Type specific carbohydrate of group A *Streptococcus*. Does not react with other Strep groups.

Form: Liquid

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: In 10 mM PBS, pH 7.2 (0.1% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -20°C.

Aliquot to avoid repeated freezing and thawing.