

Datasheet

8-Hydroxyguanosine (8OHdG) polyclonal antibody

Catalog Number: PAB13868

Regulation Status: For research use only (RUO)

Product Description: Goat polyclonal antibody raised against 8-Hydroxyguanosine (8OHdG).

Immunogen: 8-Hydroxyguanosine-conjugated.

Host: Goat

Applications: ELISA, IHC-P

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: 8-Hydroxy-2'-deoxyguanosine and 8-hydroxyguanosine.

Form: Liquid

Recommend Usage: ELISA (>1:250000)

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) (1:200)

The optimal working dilution should be determined by the end user.

Storage Buffer: In buffer containing 0.1% sodium azide

Storage Instruction: Store at 4°C. For long term storage store at -20°C.

Aliquot to avoid repeated freezing and thawing.