

Datasheet

AMY2A/AMY2B polyclonal antibody

Catalog Number: PAB13874

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against native AMY2A/AMY2B.

Immunogen: Native purified human pancreatic AMY2A/AMY2B.

Host: Rabbit

Reactivity: Human

Applications: RIA, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: Cross-reactivity against other forms of amylase has not been determined.

Form: Liquid

Recommend Usage: Western Blot (1:8000 strong reaction to 0.15ug and 0.3ug antigen)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.05% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -20°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 279|280

Gene Symbol: AMY2A

Gene Alias: AMY2, AMY2B, PA