

Datasheet

Amyloid Precursor Protein polyclonal antibody

Catalog Number: PAB13876

Regulation Status: For research use only (RUO)

Product Description: Goat polyclonal antibody raised against synthetic peptide of Amyloid Precursor Protein.

Immunogen: A synthetic peptide corresponding to amino acids 339-348 of human Amyloid Precursor Protein.

Sequence: RGR TSSKELA

Host: Goat

Reactivity: Human

Applications: ELISA, IHC-P, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: Amino acid sequence 339-348 of human amyloid precursor protein frameshift mutant.

Form: Liquid

Recommend Usage: ELISA (>1:5000)

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) (>1:300)

Western Blot (>1:3000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.1% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -20°C.

Aliquot to avoid repeated freezing and thawing.