

Datasheet

SASH1 monoclonal antibody (M02), clone X1

Catalog Number: H00023328-M02

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a partial recombinant SASH1.

Clone Name: X1

Immunogen: SASH1 (NP_056093, 1066 a.a. ~ 1175
a.a) partial recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

ENTSLQEHEGVKLGPAITRKVSCARGVDLETLTENKLH
AEGIDLTEEPYSDKHGRCGIPEALVQRYAEDLDQPER
DVAANMDQIRVKQLRKQHRMAIPSGGLTEICRKPVS

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re

(See our web site product page for detailed applications
information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG1 Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 23328

Gene Symbol: SASH1

Gene Alias: KIAA0790, RP3-323M4.1, SH3D6A,
dJ323M4, dJ323M4.1