

Datasheet

COL6 polyclonal antibody

Catalog Number: PAB14349

Regulation Status: For research use only (RUO)

Product Description: Goat polyclonal antibody raised against native COL6.

Immunogen: Native human COL6.

Host: Goat

Reactivity: Bovine, Human

Applications: Dot, ELISA, ICC, IHC-Fr, WB
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: Reacts with conformational determinants on bovine, human type VI collagen as demonstrated by ELISA. May react with type VI collagen from other species. Exhibits <10% cross reactivity with collagen type I, II, III, IV, V. This antibody has not been tested for reactivity with other ECM proteins (e.g., laminin, fibronectin).

Form: Liquid

Recommend Usage: Immunocytochemistry (1:10-1:20)
Immunohistochemistry (1:10-1:20)
The optimal working dilution should be determined by the end user.

Storage Buffer: In 100 mM BBS, pH 8.2

Storage Instruction: Store at 4°C.